



# Pre-requisites

Practice Name \_\_\_\_\_ Doctor's Signature \_\_\_\_\_

Patient's Name \_\_\_\_\_ Patient's Chart # \_\_\_\_\_

Practice Address \_\_\_\_\_

Office Phone \_\_\_\_\_ Due Date \_\_\_\_\_

## Crown and Bridge

### Full Solid Restorations

- CB 1.1 Kortz Solid Zirconia ★
- CB 1.2 Emax CAD ★
- CB 1.3 Bruxzir Solid Zirconia
- CB 1.4 Bruxzir Anterior
- CB 1.5 Nori-Take Katana ★
- CB 1.6 Lava Esthetic

### Porcelain Layered Restorations

- CB 2.1 PFM-NP Non-Precious
- CB 2.2 PFM-SP Semi-Precious 2% Au
- CB 2.3 PHN White High Noble 40% Au
- CB 2.4 PFG Yellow Gold 90% Au
- CB 2.5 PFZ Porcelain Fused to Zirconia ★

### Full Cast Metal

- CB 3.1 Full Metal Crown N.P.
- CB 3.2 Full Yellow Gold Crown 2% Au



## Removable

★ Restorations with a Green Star represent a popular choice

### Denture

- RM 1.1 Bite
- RM 1.2 Teeth Set-Up
- RM 1.3 Finish

### Flexible Partial ★

- RM 3.1 Bite
- RM 3.2 Teeth Set-Up
- RM 3.3 Finish

### Additional Services

- RM 5.1 Stayplate Flipper
- RM 5.2 Night Guard
- RM 5.3 Reline
- RM 5.4 Repair

### Framework

- RM 2.1 Frame + Bite
- RM 2.2 Teeth Set-Up
- RM 2.3 Finish

### Flexible Partial + Frame Combo

- RM 4.1 Frame + Bite
- RM 4.2 Teeth Set-Up
- RM 4.3 Finish



## Implants

**\*Please make sure to select a Crown for the Restoration from above**  
Kortz Solid Zirconia Crown is included with the price of the abutment. Any other crown will be charged at full price.

### Abutment

- IP 1.1 Titanium Abutment
- IP 1.2 Zirconia Hybrid Abutment

### System & Platform

- |  |                               |                               |                              |
|--|-------------------------------|-------------------------------|------------------------------|
| Nobel Replace<br><small>E-Series</small>             | <input type="radio"/> 3.5 NP  | <input type="radio"/> 4.3 RP  | <input type="radio"/> 5.0 WP |
| Nobel Active ★<br><small>F-Series</small>            | <input type="radio"/> 3.5 NP  | <input type="radio"/> 4.3 RP  | <input type="radio"/> 5.0 WP |
| Biomet 3i Certain<br><small>H-Series</small>         | <input type="radio"/> 3.4 NP  | <input type="radio"/> 4.1 RP  | <input type="radio"/> 5.0 WP |
| Biomet 3i OSSEOTITE<br><small>I-Series</small>       |                               | <input type="radio"/> 4.1     |                              |
| NobelBranemark<br><small>K-Series</small>            | <input type="radio"/> 3.5     | <input type="radio"/> 4.1     | <input type="radio"/> 5.1    |
| Straumann SynOcta<br><small>N-Series</small>         | <input type="radio"/> 3.5 NN  | <input type="radio"/> 4.8 RN  | <input type="radio"/> 6.5 WN |
| Zimmer Tapered Screw-Vent<br><small>R-Series</small> | <input type="radio"/> 3.5     | <input type="radio"/> 4.5     | <input type="radio"/> 5.7    |
| Astra Tech OsseoSpeed<br><small>S-Series</small>     | <input type="radio"/> 3.5/4.0 | <input type="radio"/> 4.5/5.0 |                              |
| Dentsply Friadent XiVE<br><small>T-Series</small>    | <input type="radio"/> 3.4     | <input type="radio"/> 4.5     | <input type="radio"/> 5.5    |

### Hybrid Implant Crown

- IP 1.3 Implant Crown with Xi-Base


### Screw Access Hole

- Yes
- No


### Cemented In Lab

- Yes - One Piece
- No - Two Pieces

*\*Restorations come with lab and final screw. Restorations milled from FDA Approved NT-Trading Premills. Each Trademark/Registered name is the property of their respective owner.*



## Notes



1345 N. Fitzgerald Ave. Unit F Rialto, CA 92376  
Toll Free: 888-345-0173 Phone: 909.562.0001 Web: www.HDentallab.com